

Product Datasheet

Human FH protein (orb49395)

Description	Recombinant human FH protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Ala45-Lys510
Source	E.coli
Uniprot ID	P07954
Expiration Date	6 months from date of receipt.