

Product Datasheet

Human Grancalcin protein (orb49417)

Description	Recombinant human Grancalcin protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 500mM NaCl, pH 8.5 .
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 500mM NaCl, pH 8.5 .
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Protein Sequence	Met1-Ile217
Source	E.coli
Uniprot ID	P28676
Expiration Date	6 months from date of receipt.