

Product Datasheet

Human PLA2G16 protein (orb49423)

Description	Recombinant human PLA2G16 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Asp12-Asp132
Source	E.coli
Uniprot ID	P53816
Expiration Date	6 months from date of receipt.