

Product Datasheet

Glucocorticoid Receptor beta antibody (PE) (orb494402)

Description	Glucocorticoid Receptor beta antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Glucocorticoid Receptor beta
Target	Glucocorticoid Receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04150
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.