

Product Datasheet

Human Kallikrein 3 protein (orb49490)

Description	Recombinant human Kallikrein 3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl ₂ , pH 6.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl ₂ , pH 6.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ala18-Pro261
Source	Human Cells
Uniprot ID	P07288
Expiration Date	6 months from date of receipt.