

Product Datasheet

Human KPNA2 protein (orb49497)

Description	Recombinant human KPNA2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 20% Glycerol, pH 8.0 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 1mM DTT, 20% Glycerol, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Phe529
Source	E.coli
Uniprot ID	P52292
Expiration Date	6 months from date of receipt.