


Product Datasheet

Human Legumain protein (orb49502)

Description	Recombinant human Legumain protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ile18-Tyr433
Source	Human Cells
Uniprot ID	Q99538
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Legumain
prot...