

Product Datasheet

alpha Estrogen Receptor antibody (PE) (orb495166)

Description	alpha Estrogen Receptor antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from the middle of human Estrogen Receptor alpha
Target	Estrogen Receptor alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67kDa
Uniprot ID	P03372
Expiration Date	12 months from date of receipt.