

Product Datasheet

Integrin alpha 4 antibody (PE) [Out of stock] (orb495228)

Description	Integrin alpha 4 antibody (PE) [Out of stock]
Species/Host	Rabbit
Conjugation	PE
Tested Applications	ELISA, FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Integrin alpha 4
Target	Integrin alpha 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13612
Dilution Range	WB: 1:500-2000, Flow-Cyt: 1µg/Test, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.