


Product Datasheet

Human MMP3 protein (orb49535)

Description	Recombinant human MMP3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 0.05% Brij35, 10% Glycerol, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 0.05% Brij35, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Tyr18-Cys477
Source	Human Cells
Uniprot ID	AAA36321.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of using Human MMP3 pr...