

Product Datasheet

ChRM3 antibody (PE) (orb495382)

Description	Rabbit polyclonal antibody against ChRM3 conjugated to PE
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human ChRM3
Target	Chrm3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q9ERZ3
Expiration Date	12 months from date of receipt.