

Product Datasheet

p70 S6 Kinase beta antibody (PE) (orb495461)

Description	p70 S6 Kinase beta antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from the middle of human p70 S6 Kinase
Target	RPS6KB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23443
Dilution Range	Flow-Cyt=1ug/test,IF=1:100-500
Expiration Date	12 months from date of receipt.