

Product Datasheet

GRIN2B (phospho-Tyr1472) antibody (PE) [Out of stock] (orb495903)

Description

Rabbit polyclonal antibody against

conjugated to PE

Species/Host

Rabbit

Conjugation

PE

Tested

FC, IF

Applications
Immunogen

KLH conjugated synthetic peptide derived from mouse NMDAR2B

Target

GRIN2B (phospho-Tyr1472) (PE)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A

MW

163 kDa

CAS Number

Glutamate receptor ionotropic, NMDA 2B

Uniprot ID
[Q13224](https://www.uniprot.org/entry/Q13224)
Dilution Range

Flow-Cyt: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.