

Product Datasheet

Human LpPLA2 protein (orb49591)

Description	Recombinant human LpPLA2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM HAc-NaCl, 150mM NaCl, 10% Glycerol, pH 4.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM NaAc, 150mM NaCl, 10% Glycerol, pH 5.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Phe22-Asn441
Source	Human Cells
Uniprot ID	AAH38452.1
Expiration Date	6 months from date of receipt.