

Product Datasheet

Human TEM7 protein (orb49594)

Description	Recombinant human TEM7 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Threhalose, pH 7.2.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, 5% Threhalose, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Leu19-Thr426
Source	Human Cells
Uniprot ID	AAH36059.1
Expiration Date	6 months from date of receipt.