

Product Datasheet

Acetyl-Histone H3 (K23) antibody (PE) (orb496007)

Description	Acetyl-Histone H3 (K23) antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised acetylpeptide derived from human Histone H3 around the acetylation site of K23
Target	Histone H3 (Acetyl K23)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.