

Product Datasheet

Acetyl-Histone H3 (K18) antibody (PE) (orb496008)

Description	Acetyl-Histone H3 (K18) antibody (PE)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	KLH conjugated synthesised acetylpeptide derived from human Histone H3 around the acetylation site of K18
Target	Histone H3 (Acetyl K18)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Histone H3.1
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.