

Product Datasheet

Human S100P protein (orb49620)

Description	Human S100P protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 1mM DTT, 0.1M NaCl, 20% Glycerol, pH 8.0 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 1mM DTT, 100mM NaCl, 20% Glycerol, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Lys95
Source	E.coli
Uniprot ID	P25815
Expiration Date	6 months from date of receipt.