

Product Datasheet

Human SCO1 protein (orb49621)

Description	Recombinant human SCO1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50mM PB, 1mM DTT, pH 7.2.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM PB, 1mM DTT, pH 7.2.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Protein Sequence	Gly132-Ser300
Source	E.coli
Uniprot ID	O75880
Expiration Date	6 months from date of receipt.