

Product Datasheet

Human SEMG1 protein (orb49637)

Description	Recombinant human SEMG1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Hac-NaAc, 150mM NaCl, pH 4.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Hac-NaAc, 150mM NaCl, pH 4.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Gln24-Thr402
Source	Human Cells
Uniprot ID	AAH07096.1
Expiration Date	6 months from date of receipt.