

Product Datasheet

Human NEDP1 protein (orb49639)

Description	Recombinant human NEDP1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Lys212
Source	E.coli
Uniprot ID	Q96LD8
Expiration Date	6 months from date of receipt.