

Product Datasheet

Human SPINK1 protein (orb49640)

Description	Recombinant human SPINK1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl ₂ , 1mM DTT, 0.05% Brij35, 10% Glycerol, pH 6.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM MES, 150mM NaCl, 2mM CaCl ₂ , 1mM DTT, 0.05% Brij35, 10% Glycerol, pH 6.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Asp24-Cys79
Source	Human Cells
Uniprot ID	P00995
Expiration Date	6 months from date of receipt.