

Product Datasheet

Human PPP1CC protein (orb49642)

Description	Recombinant human PPP1CC protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, pH 8.0 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 1mM DTT, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&C-His
Note	For research use only
Application notes	6His tag at the N-terminus, 6His tag at the C-terminus.
Protein Sequence	Met1-Lys323
Source	E.coli
Uniprot ID	P36873
Expiration Date	6 months from date of receipt.