

Product Datasheet

FOXL2 antibody (PE) (orb496638)

Description	Rabbit polyclonal antibody against FOXL2 conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human FOXL2
Target	FOXL2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45-50kDa
Uniprot ID	P58012
Expiration Date	12 months from date of receipt.