

Product Datasheet

Human Shh protein (orb49665)

Description	Recombinant human Shh protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 100mM NaCl, 1mM DTT, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 100mM NaCl, 1mM DTT, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Cys24-Gly197
Source	E.coli
Uniprot ID	Q15465
Expiration Date	6 months from date of receipt.