

Product Datasheet

Human THOP1 protein (orb49676)

Description	Recombinant human THOP1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB, 100mM NaCl, 10% Glycerol, pH 7.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 100mM NaCl, 10% Glycerol, pH 7.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Lys2-Cys689
Source	E.coli
Uniprot ID	P52888
Expiration Date	6 months from date of receipt.