

Product Datasheet

Human GTF2B protein (orb49686)

Description	Recombinant human GTF2B protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Protein Sequence	Met1-Leu316
Source	E.coli
Uniprot ID	Q00403
Expiration Date	6 months from date of receipt.