

## Product Datasheet

### Human PGP9.5 protein (orb49709)

|                          |                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant human PGP9.5 protein                                                                    |
| <b>Endotoxins</b>        | Less than 0.1 ng/ug (1 EU/ug)                                                                       |
| <b>Preservatives</b>     | Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 250mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.  |
| <b>Form/Appearance</b>   | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 250mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5. |
| <b>Storage</b>           | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.              |
| <b>Tag</b>               | C-His                                                                                               |
| <b>Note</b>              | For research use only                                                                               |
| <b>Application notes</b> | 6His tag at the C-terminus.                                                                         |
| <b>Protein Sequence</b>  | Met1-Ala223                                                                                         |
| <b>Source</b>            | E.coli                                                                                              |
| <b>Uniprot ID</b>        | <b>P09936</b>                                                                                       |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                      |