

Product Datasheet

Human HR6A protein (orb49713)

Description	Recombinant human HR6A protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-GST&His
Note	For research use only
Application notes	GST, 6His tag at the N-terminus.
Protein Sequence	Met1-Cys152
Source	E.coli
Uniprot ID	P49459
Expiration Date	6 months from date of receipt.