

Product Datasheet

Human GALE protein (orb49723)

Description	Recombinant human GALE protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM TrisHCl, 150mM NaCl, 2mM DTT, 1mM EDTA, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM Tris-HCl, 150mM NaCl, 2mM DTT, 1mM EDTA, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Ala348
Source	E.coli
Uniprot ID	Q14376
Expiration Date	6 months from date of receipt.