

Product Datasheet

Human ZNF100 protein (orb49749)

Description	Recombinant human ZNF100 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 50mM Imidazole, 1mM ZnCl ₂ , 30% Glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 200mM NaCl, 50mM Imidazole, 1mM ZnCl ₂ , 30% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Lys109
Source	E.coli
Uniprot ID	Q8IYN0
Expiration Date	6 months from date of receipt.