

Product Datasheet

Human betaDefensin 1 protein (orb49762)

Description	Recombinant human betaDefensin 1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 130mM NaCl, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 130mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Gly22-Lys68
Source	E.coli
Uniprot ID	P60022
Expiration Date	6 months from date of receipt.