

Product Datasheet

Mouse IL1 alpha protein (orb49778)

Description	Recombinant mouse IL1 alpha protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50mM TrisHCl, 200mM NaCl, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM Tris-HCl,200mM NaCl,pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Ser115-Ser270
Source	E.coli
Uniprot ID	P01582
Expiration Date	6 months from date of receipt.