

Product Datasheet

Complement alpha 2 Macroglobulin antibody (PE)
[Out of stock] (orb497795)

Description

Rabbit polyclonal antibody

Species/Host

Rabbit

Conjugation

PE

Immunogen

KLH conjugated synthetic peptide
derived from human
Complement C3b alpha' chain

Target

C3

Preservatives

0.01M TBS(pH7.4) with 1%
rAlbumin, 0.03% Proclin300 and
50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year.
Avoid repeated freeze/thaw
cycles. The lyophilized antibody
is stable at room temperature for
at least one month and for
greater than a year when kept at
-20°C. When reconstituted in
sterile pH 7.4 0.01M PBS or
diluent of antibody the antibody
is stable for at least two weeks at
2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

101/181kDa

Uniprot ID
P01024
Expiration Date

12 months from date of receipt.