

Product Datasheet

ERK1/2 antibody (PE) (orb499035)

Description	ERK1/2 antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human ERK1/2
Target	ERK1 + ERK2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63844
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.