

Product Datasheet

Collagen I antibody (PE) (orb499301)

Description	Collagen I antibody (PE)
Species/Host	Mouse
Reactivity	Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Collagen I
Target	Collagen I
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P02452
Dilution Range	IF=1:50-500,ICC=1:50-500
Expiration Date	12 months from date of receipt.