


Product Datasheet


LC3A antibody (orb499565)

Description	LC3A antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LC3A (1-100/121 aa)
Target	LC3A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	14/16 kDa
Uniprot ID	Q9H492
Dilution Range	WB=1:500-1000, IHC-P=1:100-500, IF=1:200-800, IHC-F=1:100-500, ICC=1:100
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Cerebrum (Rat); Lane

...



Western blot: HeLa
Cell; (Human); Primary: A...