



Product Datasheet

MEK2 antibody (orb499566)


Description	MEK2 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human MEK2 Protein (1-171/400 aa)
Target	MEK2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P36507
Dilution Range	WB=1:1000-2000, IHC-P=1:200-1000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human brain
glioma usi...



IHC-P
analysis of
rat brain
using MEK2
M...



IHC-P
analysis of
Mouse brain
using MEK2...