


Product Datasheet

ABCG2 antibody (orb499591)

Description	ABCG2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ABCG2 (561-655/655 aa)
Target	ABCG2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72 kDa
Uniprot ID	Q9UNQ0
Dilution Range	IHC-P=1:100-500, ICC=1:100, Flow-Cyt=3ug/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of A549 cell using ABCG2 po...