



Product Datasheet

AMPK Alpha-1 antibody (orb499641)


Description	AMPK Alpha-1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human AMPK alpha-1 Protein
Target	AMPK alpha-1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	64 kDa
Uniprot ID	Q13131
Dilution Range	WB=1:500-1000, IHC-P=1:200-800
Expiration Date	12 months from date of receipt.



ICC analysis of A549 cell using AMPK alp...



IHC-P analysis of Human kidney using AMP...



IHC-P analysis of rat pancreas using AMP...