



Product Datasheet

Beta-Actin antibody (orb499643)


Description	Beta-Actin antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic MAP peptide derived from human beta-Actin
Target	beta-Actin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P60709
Dilution Range	WB=1:5000-20000, IHC-P=1:200-1000, Flow-Cyt=1ug/Test, IHC-F=1:200-1000, ICC=1:100, IF=1:200-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse CO...



Immunohistochemical staining of rat colo...



Western blot analysis of 293T Cell (Huma...