

Product Datasheet

Beta Tubulin(Loading Control) antibody (orb499644)


Description

Beta Tubulin(Loading Control) antibody


Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Beta tubulin
Target	Beta tubulin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	55 kDa
Uniprot ID	P07437
Dilution Range	WB=1:5000-20000, IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:200-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of Hela cell using Beta tub...



IHC-P analysis of Mouse brain using Beta...



IHC-P analysis of rat brain using Beta t...