



Product Datasheet

Vimentin antibody (orb499662)

Description	Vimentin antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human Vimentin Protein
Target	Vimentin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	53 kDa
CAS Number	Vimentin
Uniprot ID	P08670
Dilution Range	WB=1:500-1000, IHC-P=1:100-500, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow
cytometry of
Hela.Primary
Antibody ...



IHC-P
analysis of
human
cervical
carcino...