



Product Datasheet

CK18 antibody (orb499665)


Description	CK18 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	Recombinant human CK18 protein
Target	CK18
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	48 kDa
Uniprot ID	P05783
Dilution Range	WB=1:1000-2000, IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



ICC analysis of Hela cell using CK18 monoclonal antibody



IHC-P analysis of Human colon carcinoma tissue



IHC-P analysis of Mouse liver using CK18 monoclonal antibody