



Product Datasheet

Cubilin antibody (orb4997)


Description	Cubilin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Cubilin. Please contact us for the exact immunogen sequence. The peptide is available as orb12917.
Target	Cubilin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	399 kDa
Uniprot ID	O60494
Dilution Range	IF/ICC: 1:100-600, IHC-P: 1:100-600, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



WB analysis of rat kidney (lane 1), rat ...



IF analysis of rat stomach tissue using ...



IHC-P staining of rat stomach tissue usi...