


Product Datasheet

Retinoid X Receptor Alpha antibody (orb499736)

Description	Retinoid X Receptor Alpha antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human RXR Alpha (201-300/462 aa)
Target	Retinoid X receptor alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P19793
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of Hela.Primary Antibody ...