



Product Datasheet

PARP1 antibody (orb499737)


Description	PARP1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PARP1 (581-650/1014 aa)
Target	PARP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	112 kDa
Uniprot ID	P09874
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow
cytometry of
293T.Primary
Antibody ...



IHC-P
analysis of
rat brain
using PARP1
...



Western
blot:Lane
1:K562
(Human)
Cell;La...