

Product Datasheet

NOTCH3 antibody (orb499771)

Description	NOTCH3 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NOTCH3 (2291-2321 aa)
Target	NOTCH3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	255 kDa
Uniprot ID	Q9UM47
Dilution Range	IHC-P=1:100-500, WB=1:100-500
Expiration Date	12 months from date of receipt.