

Product Datasheet

Firefly Luciferase antibody (orb499784)

Description	Firefly Luciferase antibody
Species/Host	Mouse
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Luciferase comes from the North American firefly (Firefly)
Target	Firefly Luciferase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61 kDa
Uniprot ID	P08659
Dilution Range	WB=1:500-5000
Expiration Date	12 months from date of receipt.