

Product Datasheet

H1N1 Nucleoprotein antibody (orb499850)

Description	H1N1 Nucleoprotein antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from Influenza A virus H1N1 nucleoprotein (401-498/498 aa)
Target	H1N1 nucleoprotein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55-58 kDa
Uniprot ID	P03452
Expiration Date	12 months from date of receipt.